

AFT[™] proximal humerus fracture plate

Anatomically Correct. Suture Friendly.

The **AFT[™] Proximal Humerus Fracture Plate** is intended for fractures and fracture dislocations, osteotomies, and non-unions of the proximal humerus, particularly in osteopenic bone.

- Anatomically shaped plates are provided in left and right configurations with short and long options.
- The design is optimized to enhance the surgeon's ability to employ suture.
- The plate is shaped to avoid accidental compression of the ascending branch of the anterior circumflex artery.
- The plate is undercut around the proximal periphery to allow for needle passage after the plate has been affixed to the bone.
- Diaphyseal locking screws diverge at progressively steeper angles.
- Locking screws are available in both 3.5mm and 4.5mm options.

Shoulder Options' mission is to improve patient quality of life through state-of-the-art treatments for conditions affecting the shoulder.



C. Scott Humphrey, MD
Eagle, ID

John G. Costouros, MD
Palo Alto, CA

T. Bradley Edwards, MD
Houston, TX

E. Rhett Hobgood, MD
Jackson, MS

To learn more about Shoulder Options and the AFT Proximal Humerus Fracture Plate, please visit

www.shoulderoptions.com

AFT[™] proximal humerus fracture plate

Implants

Long Fracture Plate, Left
Long Fracture Plate, Right
Short Fracture Plate, Left
Short Fracture Plate, Right
3.5mm Locking Screw, 20mm
3.5mm Locking Screw, 25mm
3.5mm Locking Screw, 30mm
3.5mm Locking Screw, 35mm
3.5mm Locking Screw, 40mm
3.5mm Locking Screw, 45mm
3.5mm Locking Screw, 50mm
3.5mm Locking Screw, 55mm
3.5mm Non-Locking Screw, 20mm
3.5mm Non-Locking Screw, 25mm
3.5mm Non-Locking Screw, 30mm
3.5mm Non-Locking Screw, 35mm
3.5mm Non-Locking Screw, 40mm
3.5mm Non-Locking Screw, 45mm
4.5mm Locking Screw, 30mm
4.5mm Locking Screw, 35mm
4.5mm Locking Screw, 40mm
4.5mm Locking Screw, 45mm
4.5mm Locking Screw, 50mm
4.5mm Locking Screw, 55mm
4.5mm Non-Locking Malleolar Screw, 30mm
4.5mm Non-Locking Malleolar Screw, 35mm
4.5mm Non-Locking Malleolar Screw, 40mm
4.5mm Non-Locking Malleolar Screw, 45mm
4.5mm Non-Locking Malleolar Screw, 50mm
4.5mm Non-Locking Malleolar Screw, 55mm

Instruments

Long Drill Guide Plate, Left
Long Drill Guide Plate, Right
Short Drill Guide Plate, Left
Short Drill Guide Plate, Right
2.7mm Drill Guide
3.0mm Drill Guide, Short
3.0mm Drill Guide, Long
3.7mm Drill Guide
2.7mm Drill
3.0mm Drill
3.7mm Drill
2.5mm Hex Driver
2.5mm Self-Retaining Hex Driver Bit
3.5mm Hex Driver
3.5mm Self-Retaining Hex Driver Bit
Proximal Depth Gauge
Distal Depth Gauge
QC/AO Ratcheting Handle
1.6x150 K-wire
Fracture Reduction Forceps
Fracture Reduction Forceps (Verbugge)

Shoulder Options, Inc.

100 E South Main St.

PO Box 1458

Waxhaw, NC 28173

Phone: (704) 512-0000

Fax: (208) 246-1466

